



COMET PATIENT MONITOR

Domain Model

Team Members:

Patrick HARING
Christian BÜRGI

Client:

Prof. Dr. Olivier BIBERSTEIN

Revision: 98


May 4, 2012

https://svn.bfh.ch/repos/projects/patmon1/trunk/doc/src/domain_model.tex?p=98

Contents

1	Domain Model	3
2	Description	4
2.1	Classes	4
2.1.1	Administrator	4
2.1.2	Doctor	4
2.1.3	Patient	4
2.1.4	Period	4
2.1.5	Device	4
2.1.6	Measurement	5
2.2	Relationships	5
2.2.1	Administrator administrates Doctor	5
2.2.2	Doctor treats Patient	5
2.2.3	Patient wears Device	5
2.2.4	Patient wears device during Period	5
2.2.5	Period has data from Device	6
2.2.6	Device is currently assigned to Period	6
2.2.7	Period has Measurement	6
2.2.8	Device takes Measurement	6

1 Domain Model



powered by Astah

2 Description

2.1 Classes

2.1.1 Administrator

The **Administrator** has the attributes **mailAddress** and **password** which are used as credentials to authenticate against the system.

2.1.2 Doctor

Each **Doctor** has an email address and a password. They are stored in the attributes **mailAddress** and **password**, which are used as credentials to authenticate against the system. **Enabled** indicates whether the account is enabled or disabled. The attributes **activated** and **activationCode** are used for the activation process. **Activated** indicates whether the **Doctor** is activated or not, **activationCode** stores the activation code which was sent to him by email. For a successful login into the system, both **enabled** and **activated** have to be true. In order to delete an account, it has to be disabled for a given time period. Therefore a timestamp called **lastStatusChange** shows when the last status change occurred. A status change is either a change of **enabled** or **activated**.

2.1.3 Patient

Each **Patient** is identified by his social security number (SSN) which is stored in the attribute **ssn**. The attributes **name**, **address** and **phone** are mandatory and will help the **Doctor** to identify the **Patient**.

2.1.4 Period

The class **Period** represents a monitoring period during which a **Device** is monitoring data on a **Patient**. The attributes **start** and **end** represent the interval of time and the attribute **frequency** represents the frequency by which the **Measurements** are taken.

2.1.5 Device

Each **Device** is identified by its serial number, which is stored in the attribute **serialNumber**.

2.1.6 Measurement

The class `Measurement` represents the `temperature` which was taken by a `Device` at a given time `timestamp`.

2.2 Relationships

2.2.1 Administrator administrates Doctor

The `Administrator` can have zero to `n` `Doctors` he administrates in the system. In the initial state of the system, there is no `Doctor`. As the system should be able to handle a big amount of `Doctors`, they are all related to the `Administrator` they were created by.

2.2.2 Doctor treats Patient

`Doctors` treat from zero to `n` `Patients`. When a `Doctor` is created, he doesn't treat any `Patient`. When doing treatments, he gets assigned all `Patients` he treats. On the other side, `Patients` can only be treated by a single `Doctor` (and there has to be a `Doctor` who treats the `Patient`).

2.2.3 Patient wears Device

When a `Device` is assigned to a `Patient`, this relationship represents that a specific `Device` is actually worn by a specific `Patient`. A `Device` can only be worn by a single `Patient` at the same time. On the other hand, a `Patient` can only wear a single `Device` at the same time. This relationship can also be represented by "`Patient` wears device during `Period`" and "`Device` is currently assigned to `Period`". Therefore it is intended that this relationship is not implemented in the actual code.

2.2.4 Patient wears device during Period

`Patients` which are monitored need to have assigned a `Period` in which the monitoring takes place. `Patients` can have multiple `Periods` if they are treated more than once, but a `Period` can only belong to a single `Patient` and the `Period` has to belong to a `Patient`.

2.2.5 Period has data from Device

In order to detect which **Device** records or has recorded the **Period**'s date, there is a relationship between **Period** and **Device**. **Periods** have exactly one **Device** since there is only one **Device** recording data for it. On the other hand, **Devices** can record data for any number of **Periods**, since they can be reconfigured to a new **Period** after they have completed monitoring for an arbitrary **Period**. The condition, that a **Device** can only be configured for a single **Period** at a time, is not directly represented in the model and has to be checked programmatically. This is done using the relationship “**Device** is currently assigned to **Period**”.

2.2.6 Device is currently assigned to Period

Periods have a second relationship to **Device**. This relationship shows, that the **Device** is configured for the **Period**. A **Device** can only be assigned to one **Period** at a time. The **Period** which the **Device** is configured to must be the same as in the relationship “**Period** has data from **Device**”.

2.2.7 Period has Measurement

Periods have a relationship to **Measurement** showing which **Measurement** has been taken for which **Period**. **Measurements** belong to exactly one **Period** but a **Period** can have any number of **Measurements** (including none).

2.2.8 Device takes Measurement

Measurements are taken by **Devices**. They belong exactly to one **Device**. On the other hand, **Devices** can take any number of **Measurements**. This relationship can be reflected using “**Period** has data from **Device**” and “**Period** has **Measurement**”. Therefore it is intended that this relationship is not implemented in the actual code.